

### **Instrument Marking**

Instruments should be embossed with

Surgins Logo

Surgins Product Code

CE Mark

### **Label Instruction**

Label should be placed on individual instrument plastic bag with

REF : Product Name and Surgins Code

LOT : Surgins PO Number which is SUR

QTY: 1

### **In final Packing,**

Individual instrument Sets should be separately packed with following label details

REF : Product Name and Surgins Code of Instrument Set

LOT : Surgins PO Number which is SUR

QTY: As per individual Sets like 6, 22 etc

### **Sharp Edges/Ends**

Please, ensure, sharp edges are protected with Plastic Tips or alternative.

**PRODUCT DESCRIPTION:**

- Re-useable Surgical Instruments, single packed in a PE Grip Seal Bag

**Batch Information:**

Lot Number: Transpack Works Order Number  
Manuf: Not Applicable.  
Expiry: Not Applicable.

**Primary Pack:**

Packed in Clear Low Density Polyethylene Grip Seal bag and sealed.  
Over-labelled according to Surgins requirements (refer label specification)

**Secondary Pack:**

Packed 10 per carton into white outer carton, CAR-0022.

**Sterilisation:**

Product is supplied non-sterile, for autoclave sterilisation prior to use.

LABEL SPECIFICATION:

8 per sheet (4 x 2 fc)

SAMPLE:



**Surgins Surgical Ltd.**

PO BOX 14530

Solihull, B91 9LH, England.

Tel.: +44(0)121 709 010

Fax: +44(0)121 705 2221

[www.surgins.co.uk](http://www.surgins.co.uk)

**REF** SUR-xx-xxxx-xx

**NON-STERILE**

**Product Description**

**LOT** Lot

**QTY** 1



Do not use if package is damaged



## FINAL PACKING

Inner Carton: CAR-0022    White outer carton, “postal” style box.

### Packing Instructions:

- Place a layer of bubble wrap in the base of the box.
- Pack the product in the box, alternating the handle end of the product in the box.  
Quantity 10 per box.
- Pack 1 x Instruction / Care leaflet in the box on top of the instruments.
- Fill the carton with bubble wrap to ensure the product enclosed is protected and does not move excessively.
- Apply 1 x label (see below) to the carton across the opening flap of the carton, with batch details visible on the base. This label will perform the function of a tamper-evident label.

Box Quantity: 10

Labels:            1 Label only to be applied to the top long face in a central position.

Label Size:    A4 8 per sheet.